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
EN

COMPOSITION

Hyaluronic Acid gel 18 mg
Phosphate buffer pH 7.2 q.s.
One syringe contains 0.6 ml of **Juvéderm™ REFINE™**.

DESCRIPTION

Juvéderm™ REFINE™ is a sterile pyrogen-free physiological solution of cross-linked hyaluronic acid which is not of animal origin. The gel is presented in a graduated pre-filled disposable syringe. Each box contains two syringes of **Juvéderm™ REFINE™**, 4 30G1/2 disposable needles reserved for injection of **Juvéderm™ REFINE™**, an instruction leaflet and a set of labels in order to ensure traceability.

STERILISATION

The content of the **Juvéderm™ REFINE™** syringes is sterilised by moist heat.

The 30G1/2 needles are sterilised by radiation.

INDICATIONS

Juvéderm™ REFINE™ is an injectable implant indicated for filling superficial skin depressions by injection in the superficial dermis.

CONTRA-INDICATIONS

- Do not inject **Juvéderm™ REFINE™** into the eyelids. The application of **Juvéderm™ REFINE™** in the bags under the eyes is reserved to specialists specifically trained in this technique and having a sound knowledge of the physiology for this particular area.
- Do not inject into the blood vessels (intravascular).
- Do not overcorrect.

- **Juvéderm™ REFINE™** must not be used in :
 - patients who tend to develop hypertrophic scarring
 - patients who are known to be hypersensitive to hyaluronic acid.
 - women who are pregnant or breastfeeding.
 - children.

- **Juvéderm™ REFINE™** must not be used in areas presenting cutaneous inflammatory and/or infectious processes (acne, herpes ...).
- **Juvéderm™ REFINE™** must not be used in association with laser therapy, chemical peeling or dermabrasion. For surface peels, it is recommended not to inject the product if the inflammatory reaction generated is significant.

PRECAUTIONS FOR USE

- **Juvéderm™ REFINE™** is not indicated for injections other than intra-dermal injections.
- As a matter of general principle, injection of a medical device is associated with a risk of infection.
- There are no available clinical data (efficiency, tolerance) about injection of **Juvéderm™ REFINE™** into an area which has already been treated with another filling product. It is recommended not to inject it in site which has been treated with a permanent implant.
- No clinical data is available regarding the efficiency and tolerance of **Juvéderm™ REFINE™** injections in patients having a history of or currently suffering from autoimmune disease. Depending on the nature of the disease, the medical practitioner shall therefore decide on the administration and associated treatment on an individual case-by-case basis. He/ she shall also ensure the specific monitoring of this type of patients. In particular, it is recommended that these patients undergo a preliminary dual test, and to refrain from injecting the product if the disease is active.
- There is no available clinical data concerning the tolerance of the **Juvéderm™ REFINE™** injection in patients presenting a history of severe multiple allergies or anaphylactic shock. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the allergy, and shall also ensure the specific monitoring of these at-risk patients. In particular, the decision may be taken to propose a double test or suitable preventive treatment prior to any injection.
- Patients showing a history of streptococcal disease (recurrent sore throat, acute rheumatic fever) shall be subjected to a dual test before any injection is given to them. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.
- Patients on anti-coagulation medication (anticoagulants, aspirin, or nonsteroidal antiinflammatory drugs) must be warned of the potential increase risks of haematomas and bleeding during injection.
- There is no data available regarding the safety of injecting greater amount than 20 mL of **Juvéderm™ REFINE™** per 60 kg (130 lbs) body mass per year.
- Patients should be recommended not to apply make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight, UV light, freezing temperatures or using saunas or Turkish baths for the two weeks after the injection.
- If the 30G1/2 needle is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the needle.
- The composition of this product is compatible with fields used for magnetic resonance imaging.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. **Juvéderm™ REFINE™** should never therefore be placed in contact with these substances or with medical-surgical instrumentation which has been treated with this type of substance.

UNDESIRABLE EFFECTS

Medical practitioners must inform the patient that there are potential side effects associated with implantation of this device, which may occur immediately or may be delayed. These include (non-exhaustive list) :

- Inflammatory reactions (redness, oedema, erythema, ...) which may be associated with itching, pain on pressure, occurring after the injection. These reactions may last for a week.
- Haematomas.
- Induration or nodules at the injection site.
- Staining or discolouration of the injection site.
- Poor effect or weak filling effect.
- Cases of necroses in the glabellar region, abscesses, granuloma, and immediate or delayed hypersensitivity after hyaluronic acid injections have been reported. It is therefore advisable to take these potential risks into account.

Patients must report inflammatory reactions which persist for more than one week or any other secondary effect which develops, to their medical practitioner as soon as possible. The medical practitioner should treat these as appropriate.

Any other undesirable side effects associated with injection of **Juvéderm™ REFINE™** must be reported to the distributor and/ or to the manufacturer.

METHOD OF USE – POSOLOGY

This product is designed to be injected into the dermis by a medical practitioner. As precision is essential to a successful treatment, the product must be used by medical practitioners who have undertaken specific training in injection techniques for filling.

- **Juvéderm™ REFINE™** is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity and performance of the product and it can therefore no longer be assured.
- Before starting treatment patients should be informed of the product's indications, contra-indications, incompatibilities and potential undesirable effects.
- The area to be treated should be disinfected thoroughly prior to the injection.

Remove tip cap by pulling it straight off the syringe as shown in (fig. 1). Then firmly push the needle provided in the box (fig. 2) into the syringe, screwing it gently clockwise. Twist once more until it is fully locked and has the needle cap in the position shown in fig. 3. If the needle cap is positioned as shown in fig. 4, it is incorrectly attached.

Next, remove the protective cap by holding the body of the syringe in one hand, the protective cap in the other, as shown in fig. 5, and pulling the two hands in opposite directions. Inject slowly.

Failure to comply with these precautions could cause a disengagement of the needle and/or product leakage at luer-lock level.

- The amount injected will depend on the areas which are to be corrected.
- It is important to massage the area treated after the injection in order to ensure that the substance has been uniformly distributed.

WARNINGS

- Confirm the expiry date on the product label.
- Do not re-use. Sterility of this device cannot be guaranteed if the device is re-used.
- Do not re-sterilise.
- For the needles (C E 0123 TSK Laboratory, Japan):

EC-Representative : Emergo Europe
Molenstraat 15
2513 BH The Hague (NL)

- Used needles must be thrown away in the appropriate containers. Do the same for the syringe. Please consult the current applicable directives to ensure their correct elimination.
- Never try to straighten a bent needle; throw it away and replace it.

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