Juvéderm[®]
VOLBELLA[®]
with Lidocaine

2 x 1mL

 **ALLERGAN**



72525ZQ11
Revision 2015-08-26

ALLERGAN

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COMPOSITION

Hyaluronic acid gel	15 mg
Lidocaine hydrochloride	3 mg
Phosphate buffer pH 7.2 q.s.	1 mL

One syringe contains 1mL of **Juvéderm® VOLBELLA® with Lidocaine**.

DESCRIPTION

Juvéderm® VOLBELLA® with Lidocaine is a sterile pyrogen-free physiological solution of cross-linked hyaluronic acid gel which is not of animal origin. The gel is presented in a graduated pre-filled disposable syringe. Each box contains two 1mL **Juvéderm® VOLBELLA® with Lidocaine** syringes, 4 single-use 30G1/2" sterile needles specifically intended for injecting **Juvéderm® VOLBELLA® with Lidocaine**, an instruction leaflet and a set of labels in order to ensure traceability.

STERILISATION

The contents of the **Juvéderm® VOLBELLA® with Lidocaine** syringes is sterilised by moist heat. The 30G1/2" needles are sterilised by radiation.

INDICATIONS

- **Juvéderm® VOLBELLA® with Lidocaine** is an injectable gel implant intended for enhancement and pouting of the lips to correct volume loss and to treat perioral lines and oral commissures.
- **Juvéderm® VOLBELLA® with Lidocaine** is intended to be used via superficial or mid-dermis injection or lips mucosa by an authorized medical practitioner.
- The presence of lidocaine is meant to reduce the patient's pain during treatment.

CONTRA-INDICATIONS

- Do not inject **Juvéderm® VOLBELLA® with Lidocaine** into the eyelids.
- Do not inject into the blood vessels (intravascular).
- Do not overcorrect.
- **Juvéderm® VOLBELLA® with Lidocaine** must not be used in :
 - patients suffering from untreated epilepsy;
 - patients who tend to develop hypertrophic scarring;
 - patients who are known to be hypersensitive to hyaluronic acid;
 - patients with known hypersensitivity to lidocaine or to amide-type local anaesthetics;
 - patients suffering from porphyria;

- women who are pregnant or breastfeeding;
- children.

• **Juvéderm® VOLBELLA® with Lidocaine** must not be used in areas presenting cutaneous inflammatory and/or infectious processes (acne, herpes...).

• **Juvéderm® VOLBELLA® with Lidocaine** should not be used simultaneously with laser treatment, deep chemical peels or dermabrasion. For surface peels, it is recommended not to inject the product if the inflammatory reaction generated is significant.

PRECAUTIONS FOR USE

- **Juvéderm® VOLBELLA® with Lidocaine** is indicated only for injections for intra-dermal injections and injections in mucous membrane of the lips.
- The maximum volume that may be injected for initial and top-up treatments combined is 2 mL for lips injection, 1.5 mL for oral commissure injection and 1.2 mL for perioral line injection.
- The safety and efficacy of **Juvéderm® VOLBELLA® with Lidocaine** for enhancement and pouting of the lips has not been established in patients under the age of 21 years.
- The safety and effectiveness for the treatment of other anatomic areas, such as periorbital area (under-eye area and crow's feet) or glabella region, have not been established in controlled clinical studies.
- As a matter of general principle, injection of a medical device is associated with a risk of infection.
- There is no available clinical data (efficiency, tolerance) about injection of **Juvéderm® VOLBELLA® with Lidocaine** into an area which has already been treated with another filling product. It is recommended not to inject it in site which has been treated with a permanent implant.
- No clinical data is available regarding the efficiency and tolerance of **Juvéderm® VOLBELLA® with Lidocaine** injections in patients having a history of or currently suffering from autoimmune disease. Depending on the nature of the disease, the medical practitioner shall therefore decide on the administration and associated treatment on an individual case-by-case basis. He/she shall also ensure the specific monitoring of this type of patients. In particular, it is recommended that these patients undergo a preliminary dual test, and to refrain from injecting the product if the disease is active.
- There is no available clinical data concerning the tolerance of the **Juvéderm® VOLBELLA® with Lidocaine** injection in patients presenting a history of severe multiple allergies or anaphylactic shock. The medical practitioner must therefore decide on the indication

according to the individual case, depending on the nature of the allergy, and must ensure that there is individual surveillance of these patients who are at risk. In particular, the decision may be taken to propose a double test or preventive adapted treatment previously to any injection.

- Patients showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a dual test before any injection is given to them. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.

- Patients on anti-coagulation medication (anticoagulants, acetylsalicylic acid, or nonsteroidal anti-inflammatory drugs) must be warned of the potential increased risks of haematomas and bleeding during injection.

- There is no data available regarding the safety of injecting greater amount than 20 mL of **Juvéderm® VOLBELLA® with Lidocaine** per 60 kg (130 lbs) body mass per year.

- The combination of **Juvéderm® VOLBELLA® with Lidocaine** with certain drugs that reduce or inhibit hepatic metabolism (cimetidine, beta-blockers, etc.) is inadvisable.

- **Juvéderm® VOLBELLA® with Lidocaine** should be used with caution in patients showing symptoms of cardiac conduction disorders.

- The safety in patients with known susceptibility to keloid formation, and pigmentation disorders has not been studied.

- Please recommend the patient not to use any make up during the 12 hours that follow the injection treatment and to avoid any extended exposure to the sun, UV rays and temperatures below 0°C, as well as any sauna or hammam session during the two weeks that follow the injection treatment.

- If the needle is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the needle.

- Athletes should be made aware that this product contains an active principle that may produce a positive result in anti-doping tests.

- Medical practitioner must take into account the fact that this product contains lidocaine.

- The composition of this product is compatible with fields used for magnetic resonance imaging.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. **Juvéderm® VOLBELLA® with Lidocaine** should never therefore be placed in contact with these substances or with medical-surgical instrumentation which has been treated with this type of substance.

UNDESIRABLE EFFECTS

The patients must be informed that there are potential side effects associated with implantation of this device, which may occur immediately or may be delayed. These include (non-exhaustive list) :

- Inflammatory reactions (redness, oedema, erythema, ...) which may be associated with itching, pain on pressure, occurring after the injection. These reactions may last for a week. In particular, it has to be noticed that injection in the mucous membrane may cause more oedema and bruising due to the specific physiology of these tissues. Besides, a preventive anti-inflammatory treatment by a medical practitioner can be recommended.

- Haematomas.

- Induration or nodules at the injection site.

- Staining or discolouration of the injection site.

- Poor effect or weak filling effect.

- Cases of necrosis, abscesses, granuloma and immediate or delayed hypersensitivity after hyaluronic acid and/or lidocaine injections have been reported. It is therefore advisable to take these potential risks into account.

- Patients must report inflammatory reactions which persist for more than one week, or any other side effect which develops, to their medical practitioner as soon as possible. The medical practitioner should use an appropriate treatment.

- Any other undesirable side effects associated with injection of **Juvéderm® VOLBELLA® with Lidocaine** must be reported to the distributor and/or to the manufacturer.

METHOD OF USE – POSOLOGY

- This product is designed to be injected into the dermis or the mucous membrane of the lips by an authorized medical practitioner in accordance with local applicable regulation.

- Use of the supplied 30G 1/2" needle is recommended. However, depending on the medical practitioner's preferred injection technique, it is possible to use a 30G sterile cannula (please refer to the list hereunder). Choice of cannula length is determined by the user according to his/her injection technique.

Material Number	Description
94323/ HPC30019ACSH	Easyflow System-20* cannula 30G x 19mm.
94324/ HPC30025ACSH	Easyflow System-20* cannula 30G x 25mm.

• Contra-Indications, Method of use, Precautions for use and Warnings defined for the needle in this leaflet apply also to the cannula referenced above if used with this product.

The technicality of the latter being essential to the success of the treatment, the device must be used by medical practitioners having undertaken specific training in injection techniques for filling skin depressions and for enhancement and pouting of the lips.

• **Juvéderm® VOLBELLA® with Lidocaine** is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity and performance of the product and it can therefore no longer be assured.

• Before starting treatment patients should be informed of the indications for the device, its contra-indications, incompatibilities and potential undesirable effects.

• The area to be treated should be disinfected thoroughly prior to the injection.

• Remove tip cap by pulling it straight off the syringe as shown in fig. 1. Then firmly push the needle provided in the box (fig. 2) into the syringe, screwing it gently clockwise. Twist once more until it is fully locked and has the needle cap in the position shown in fig. 3. If the needle cap is positioned as shown in fig. 4, it is incorrectly attached.

Next, remove the protective cap by holding the body of the syringe in one hand, the protective cap in the other, as shown in fig. 5, and pulling the two hands in opposite directions.

Inject slowly.

Failure to comply with these precautions could cause a disengagement of the needle and/or product leakage at luer-lock level.

• The amount injected will depend on the areas which are to be corrected. In particular, injection of an excessive volume can be at the origin of severe oedema.

• It is important to massage the area treated after the injection in order to ensure that the substance has been uniformly distributed.

WARNINGS

• Check the expiry date on the product label.
• Do not re-use. Sterility of this device cannot be guaranteed if the device is re-used.

• Do not re-sterilise.

• For the needles (**CE** 0123 TSK Laboratory, Japan):

EC-Representative : Emergo Europe

Molenstraat 15

2513 BH The Hague (NL)

- Used needles must be thrown away in the appropriate containers. Do the same for the syringes. Please consult the current directives in force to ensure their correct elimination.

- Never try to straighten a bent needle, throw it away and replace it.

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