27

73

70

68

1 x 1mL

72448JR12
Revision 2016-05-03

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




EN

COMPOSITION

Hyaluronic acid	13.5 mg
Mannitol	9 mg
Phosphate buffer pH 7.2 q.s.	1 mL

One syringe contains 1 mL of **Juvéderm® HYDRATE**.

DESCRIPTION

Juvéderm® HYDRATE is a sterile pyrogen-free physiological solution of hyaluronic acid which is not of animal origin containing 0.9% mannitol. The gel is presented in a graduated, pre-filled, disposable syringe. Each box contains one 1mL **Juvéderm® HYDRATE** syringe, 2 sterile 30G1/6" needles to be used only for injecting **Juvéderm® HYDRATE**, an instruction leaflet and a set of labels in order to ensure traceability.

STERILISATION

The contents of the **Juvéderm® HYDRATE** syringes is sterilised by moist heat.

The 30G1/6" needles are sterilised by ethylene oxide.

INDICATIONS

Juvéderm® HYDRATE is a solution indicated for use in improving skin hydration and elasticity by multi-injection into the dermal epidermal junction and into the superficial dermis.

CONTRA-INDICATIONS

- Do not inject **Juvéderm® HYDRATE** into the eyelids. The application of **Juvéderm® HYDRATE** in the bags under the eyes is reserved to specialists specifically trained in this technique and having a sound knowledge of the physiology for this particular area.
- Do not inject into the blood vessels (intravascular). Intravascular injection may lead to embolization, occlusion of the vessels, ischemia or infarction.
- **Juvéderm® HYDRATE** must not be used in:
 - Patients who tend to develop hypertrophic scarring;
 - Patients with known hypersensitivity to hyaluronic acid and/or to gram positive bacterial proteins as hyaluronic acid is produced by Streptococcus type bacteria or mannitol;
 - Women who are pregnant or breastfeeding;
 - Children.
- **Juvéderm® HYDRATE** must not be used in areas presenting cutaneous inflammatory and/or infectious processes (acne, herpes, etc.).
- **Juvéderm® HYDRATE** should not be used simultaneously with laser treatment, deep chemical peels or dermabrasion. For surface peels, it is recommended not to inject the product if the inflammatory reaction generated is significant.

PRECAUTIONS FOR USE

- **Juvéderm® HYDRATE** is indicated only for intradermal injections.
- As a matter of general principle, injection of a medical device is associated with a risk of infection. Standard precautions associated with injectable materials shall be followed.

- There is no available clinical data about injection of **Juvéderm® HYDRATE** into an area which has already been treated with a non-ALLERGAN dermal filler.
- It is recommended not to inject into a site which has been treated with a permanent implant.
- No clinical data is available regarding the efficiency and tolerance of **Juvéderm® HYDRATE** injections in patients having a history of, or currently suffering from, autoimmune disease or autoimmune deficiency or being under immunosuppressive therapy. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the disease and its corresponding treatment, and shall also ensure the specific monitoring of these patients. In particular, it is recommended that these patients undergo a preliminary skin testing for hypersensitivity, and to refrain from injecting the product if the disease is active.
- There is no available clinical data concerning the tolerance of **Juvéderm® HYDRATE** injection in patients presenting a history of severe and/or multiple allergies. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the allergy, and shall also ensure the specific monitoring of these at-risk patients. In particular, the decision may be taken to propose a skin testing for hypersensitivity or suitable preventive treatment prior to any injection. In case of history of anaphylactic shock, it is recommended not to inject the product.
- Patients showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a skin testing for hypersensitivity before any injection is administered. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.
- Patients on anti-coagulation medication or using substances that can prolong bleeding (warfarin, acetylsalicylic acid, nonsteroidal anti-inflammatory drugs, or other substances known to increase coagulation time such as herbal supplements with garlic or ginkgo biloba, etc.) must be warned of the potential increased risks of bleeding and haematomas during injection.
- Please recommend that the patient not use any makeup during the 12 hours following the injection treatment and that any extended exposure to the sun, UV rays and temperatures below 0°C be avoided, as well as any sauna or hammam sessions during the two weeks following the injection treatment.
- The composition of this product is compatible with fields used for magnetic resonance imaging.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. **Juvéderm® HYDRATE** should therefore never be placed in contact with these substances or with medical-surgical instrumentation which has been treated with this type of substance.

UNDESIRABLE EFFECTS

The patients must be informed that there are potential side effects associated with implantation of this product, which may occur immediately or may be delayed. These include, but are not limited to:

- Inflammatory reactions (redness, oedema, erythema, etc.) which may be associated with itching and/or pain on pressure

and/or paresthesia, occurring after the injection. These reactions may last for a week.

- Haematomas.
- Induration or nodules at the injection site.
- Staining or discolouration of the injection site might be observed, especially when hyaluronic acid based gel is injected too superficially and/or in thin skin (Tyndall effect).
- Poor effect or weak treatment effect.
- Rare but serious adverse events associated with intravascular injection of hyaluronic acid in the face and tissue compression have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis and damage to underlying structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in the vision, signs of stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate medical practitioner specialist should an intravascular injection occur. Abscesses, granuloma and immediate or delayed hypersensitivity after hyaluronic acid injections have also been reported. It is therefore advisable to take these potential risks into account.
- Patients must report inflammatory reactions which persist for more than one week, or any other side effect which develops, to their medical practitioner as soon as possible. The medical practitioner should use an appropriate treatment.
- Any other undesirable side effects associated with injection of **Juvéderm® HYDRATE** must be reported to the distributor and/or to the manufacturer.

METHOD OF USE - POSOLOGY

- This product is designed to be injected by an authorized medical practitioner in accordance with local applicable regulation. In order to minimize the risks of potential complications and as precision is essential to a successful treatment, the product should be only used by medical practitioners who have appropriate training, experience and who are knowledgeable about the anatomy at and around the site of injection.
- **Juvéderm® HYDRATE** is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity and performance of the product and it can therefore no longer be assured.
- Prior to treatment, medical practitioners shall inform their patients about the product's indications, contra-indications, incompatibilities and potential associated undesirable effects/risks and ensure that patients are aware of signs and symptoms of potential complications.
- If storage in a refrigerator (see storage conditions) leave the product to come to room temperature before performing the injection.
- The area to be treated should be disinfected thoroughly prior to the injection. In particular, for large surface treatment, the whole area shall be considered (e.g. full face).
- Remove tip cap by pulling it straight off the syringe as shown in fig. 1. Hold the syringe body, firmly insert the needle provided in the package (fig. 2). Firmly attach the needle turning it gently clockwise as shown fig. 2 until you get it well engaged into the syringe Luer Lock system. Check the needle visually according to

figs. 3 and 4. Holding the syringe body in one hand and needle cap in the other, pull the two hands in opposite direction to remove it, as shown in fig. 5.

Prior to injecting, depress the plunger rod until the product flows out of the needle.

Inject slowly and apply the least amount of pressure necessary. If the needle is blocked, do not increase the pressure on the plunger rod. Instead, stop the injection and replace the needle. Failure to comply with these precautions could cause a disengagement of the needle and/or product leakage at luer-lock level and/or increase the risk of vascular compromise.

- If immediate blanching occurs at any time during the injection, the injection should be stopped and appropriate action taken such as massaging the area until its return to a normal color.
- The amount injected will depend on the areas to be treated based on the experience and the injection technique of the medical practitioner. Injection of an excessive volume can be at the origin of some side effects such as tissue necrosis and oedema.
- It is important to massage the area treated after the injection in order to ensure that the substance has been uniformly distributed.
- It is recommended not to repeat treatments more frequently than every 15 days.

WARNINGS

- Check the expiry date on the product label.
- In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe.
- Do not re-use. Sterility of this device cannot be guaranteed if the device is re-used.
- Do not re-sterilise.
- For the needles (**CE0476 GALLINI Italy**):
 - Used needles must be thrown away in the appropriate containers. Do the same for the syringes. Please consult the current applicable directives to ensure their correct elimination.
 - Never try to straighten a bent needle, throw it away and replace it.

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