Strattice[™] Reconstructive Tissue Matrix

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Device Description

Strattice™ Reconstructive Tissue Matrix (Strattice TM or Strattice surgical mesh) is a surgical mesh that is derived from porcine skin and is processed and preserved in a patented aqueous phosphate buffered solution containing matrix stabilizers. Strattice TM is intended to perform as a surgical mesh for soft tissue repair while presenting a scaffold to the patient's tissue. The structural properties minimize tissue attachment to the mesh. The Strattice surgical mesh consists of a sterilized sheet of processed porcine dermis provided in prescribed different sizes, dimensions, and thicknesses and packaged in a double pouch.

Use of the surgical mesh provides for a strong and biocompatible implant and will incorporate into the patient's tissue with associated cell and microvascular ingrowth.

Animal studies show a low incidence of adhesion to the Strattice surgical mesh based on observation of minimal visceral tissue attachment.

Indications

Strattice surgical mesh is intended for use as a soft tissue patch to reinforce soft tissue where weakness exists and for the surgical repair of damaged or ruptured soft tissue membranes. Indications for use include the repair of hernias and/or body wall defects which require the use of reinforcing or bridging material to obtain the desired surgical outcome.

It is indicated to be used to reconstruct, to recontour and to reform the host's human soft connective tissue particularly where loss of tissue has occurred and as a supporting tissue in surgical procedures such as abdominal wall hernia repairs and breast reconstruction.

Strattice surgical mesh is supplied sterile and is intended for single patient one-time use.

Contraindications

- This surgical mesh is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.
- Polysorbate 20 is a component of the aqueous phosphate buffered solution and therefore Strattice surgical mesh should not be used in patients with a known sensitivity to this material.

Warnings

- **Do not resterilize.**
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised. Do not use if the temperature monitoring device does not display "OK".
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practices and applicable national and regional environmental laws on disposal of packaging and biological materials.
- Strattice surgical mesh cannot be reused once it has been removed from the packaging and/or is in contact with a patient without increased risk of patient-to-patient contamination and subsequent infection.

Precautions

- Discard surgical mesh if handling has caused possible damage or contamination.
- Discard surgical mesh if it is past the use-by-date of the product (indicated as 4 digit year, 2 digit month, and 2 digit day [YYYY-MM-DD]).
- Ensure that the surgical mesh is put into a sterile basin and covered with room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes prior to implantation.
- Place surgical mesh in maximum possible contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.
- The Strattice surgical mesh should be hydrated and moist when the package is opened. If the Strattice surgical mesh is dry, do not use.
- In significantly contaminated or infected cases, utilize bioburden-reducing techniques to minimize contamination levels at the surgical site, including, but not limited to, appropriate drainage, debridement, negative pressure therapy, and/or antimicrobial therapy prior and in addition to implantation of the Strattice surgical mesh. (Presence of a significant microbial load may impact overall performance of the surgical mesh.)
- Large hernia defects and a bridging mesh technique are risk factors for hernia recurrence. Likewise, in large abdominal wall defect cases where midline fascial closure cannot be obtained, with or without separation of components techniques, utilization of Strattice TM in a bridged fashion is associated with a higher risk of hernia recurrence than when used to reinforce fascial closure.
- Certain considerations should be made in order to reduce the risk of adverse events when performing surgical procedures using a surgical mesh such as Strattice TM. Please see the following sections for more information: Product/Patient Selection, Technique Guidance, and Post-Operative Care.

Adverse Events

Potential adverse events are those typically associated with surgical mesh materials and/or their implantation procedures including, but not limited to, infection, foreign body response, hematoma, seroma formation, failure to integrate, recurrence of tissue defect, bulging, fistula formation, lack of tissue perfusion, inflammation, wound dehiscence and adhesion formation.

Storage

- The Strattice surgical mesh is a sterile medical device that should be stored in a clean, dry location at -8°C to 30°C.
- Refer to the temperature monitor located on the product carton to ensure that product has been stored within its temperature limits. Only use the product if the included temperature monitor displays "OK" on the screen. If screen displays anything other than "OK," do not use the product*.
- It is to be stored in its original packaging.
- The use-by-date of the product is indicated as 4 digit year, 2 digit month, and 2 digit day (YYYY-MM-DD).

Sterilization

This product is sterilized by E-Beam.

Product/Patient Selection

The following considerations are intended to serve only as general guidelines. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using Strattice surgical mesh.

- Carefully consider the risk/benefit balance of using a surgical mesh such as Strattice™ in patients with significant co-morbidities; including but not limited to, obesity, smoking, diabetes, immunosuppression, malnourishment, poor tissue oxygenation (such as COPD), and pre- or post-operative radiation.
- For breast reconstruction procedures, carefully consider patient factors (e.g. skin flap size and viability) and product sizing when selecting a surgical mesh such as Strattice™. It is important to create space for the tissue expander/breast implant that precisely fits the tissue expander/breast implant. This will help to avoid any potential space for fluid accumulation within the implant pocket. With proper product sizing the skin flaps should precisely fit over the muscle/Strattice™ layer without redundancy and without excessive tension.

Instructions For Preparing Strattice™ For Surgical Use

These instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using Strattice surgical mesh.

Required Materials

- Sterile forceps.
- Soaking fluid: room temperature sterile saline or room temperature sterile lactated Ringer's solution.
- One sterile basin per piece of the Strattice surgical mesh.

Preparation

1. Open the carton and remove the pouch.
2. Peel open the outer pouch and remove the inner pouch using aseptic technique. The inner pouch is sterile and may be placed directly into the sterile field.
3. Open the inner pouch carefully and aseptically remove the surgical mesh. Always use sterile gloves or forceps when handling Strattice surgical mesh.
4. Place the surgical mesh in a sterile basin and submerge in room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes.
5. Keep the surgical mesh in the room temperature sterile solution until ready for

implantation. Surgical mesh can be retained in sterile solution for a maximum of 4 hours.

Implantation

1. Prepare the surgical site using standard techniques.
2. Strattice surgical mesh may be folded, trimmed or cut as required to fit the surgical site using aseptic technique, ensuring allowance for overlap.
3. Transfer Strattice surgical mesh to the surgical site using sterile gloves or forceps.
4. Suture Strattice surgical mesh into place.
5. Complete the standard surgical procedure.
6. Discard any unused portions of Strattice surgical mesh as per institutional procedures.

The following suggestions are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using Strattice surgical mesh.

Technique Guidance — Breast Reconstruction

- As with any surgical implant, careful aseptic technique should be practiced and contact of the mesh with the patient's skin should be minimized.
- The quality and adequate thickness of the skin flap should be carefully assessed to determine the best method of breast reconstruction for the patient. Carefully assess the skin flaps to ensure that they are well perfused prior to proceeding with breast reconstruction using Strattice TM.
- Strattice TM should be covered with healthy, viable tissue; portions of the skin flap that appear ischemic or necrotic should be addressed using clinical judgment.
- Care should be taken to avoid excessive tension on the skin flaps at the time of closure.
- The subcutaneous space created for the tissue expander/breast implant should precisely fit the tissue expander/breast implant so as to avoid any potential space for fluid accumulation within the implant pocket. The skin flaps should precisely fit over the muscle/Strattice TM layer without redundancy and without excessive tension. Additional caution is recommended for larger implants (>500 ml) as the weight and size of the implant may contribute to the tension on the skin closure.
- Drains should be placed in the tissue expander/breast implant space as well as in the space between the skin flap and Strattice TM so as to reduce fluid accumulation.

Technique Guidance — Abdominal Wall Repair

- As with any surgical implant, careful aseptic technique should be practiced and contact of the mesh with the patient's skin should be minimized.
- Tension and suture placement are application dependent. For hernia repair



applications, surgical experience with biological surgical mesh indicates that suturing the Strattice™ under physiologic tension with a minimum of 3cm–5cm overlap, or as much as required to reach healthy adjacent tissues, may produce improved outcomes.

- Use of permanent sutures is recommended.

General Post-Operative Care

- Proper post-operative care may include the use of appropriate drainage and negative pressure therapy.
- Drains should be left in place until the amount of drainage reaches a level less than 30 ml per 24 hr period for 3 consecutive days. Typical timeline for drain removal is generally 1 to 2 weeks, but should be customized to the clinical circumstance.
- As with any post-operative care, aseptic technique should be practiced when required to minimize contamination to the surgical wound.

Definitions

 PLIABLE	 FIRM	 OPEN HERE
 CONTENTS	 INNER POUCH	 OUTER POUCH
 DO NOT USE IF PACKAGE IS DAMAGED	 CORRECT SETTING	 *INCORRECT SETTING
 X-SHAPED INCISION FOR STOMA REINFORCEMENT	 STERILIZED USING IRRADIATION	 FOR SINGLE USE ONLY
 DO NOT RE-STERILIZE	 ATTENTION, SEE INSTRUCTIONS FOR USE	 CATALOG NUMBER
 USE-BY DATE	 TEMPERATURE LIMITATIONS	 BATCH CODE

Customer Support

*Call customer support on day of receipt for immediate return authorization.

For product complaints and potential adverse events, please contact your local Sales Representative or Customer Support.

This product and certain methods are covered by patents and patents pending, including: US 8,905,826; US 8,735,054; US 8,469,779; US 8,323,352; US 8,007,531; and US 7,476,249.

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© LifeCell Corporation 2016
LifeCell Corporation
KCI Medical Products, Ltd.
KCI Medical Australia PTY, Ltd.